

## ELECTRONIC BODY FAT SCALE



**MODEL: R-GL-F8**

## INSTRUCTION SHEET

### PRODUCT INFORMATION:

It is important information for us to learn about our health situation through body fat and water percentage. This product uses BIA (Bioelectric Impedance Analysis) technology to determine through body fat and water percentage. When you stand barefoot on metal foot plates, some weak electric be collected as data (using Multi-frequency Impedance Measuring and Bars-Capacity Sensory Technology) and processor in this scale will use Population Body Composition Model Equations to determine your body fat, water, muscles and bones percentage with the collected electric data together with your input as age, height and gender

### REFERENCE TABLE:

**1.BF & BW comparison sheet**

